x-epil[®]

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x-epil[®]

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
INTENDED USE (XE9406)

The test is a pregnancy test. It is used for the qualitative detection of hCG in human urine to aid in the diagnosis of pregnancy. The test is for self-testing use. It is automated. For in vitro diagnostic use only.
For self-testing in vitro diagnostic use only. Do not swallow any part of the test. Not to be taken internally.

MATERIALS:

Materials Provided

One digital rapid pregnancy test in foil pouch.

Materials Required But Not Provided

- Timer
- Urine Collection Container

PRECAUTIONS

- For self-testing in vitro diagnostic use only. Do not swallow any part of the test. Not to be taken internally.
- Use the test only once. Do not dismantle and touch the absorbent tip of the test.
- Only use and store in conditions listed below:
 - Storage conditions: 2–30°C (35.6–86°F), _ 90% RH, 86–106 kPa;
 - Operating conditions: 15–30°C (59–86°F), _ 90% RH, 86–106 kPa.
- This product incorporates batteries, used materials should be discarded according to local regulations.
- Pollution degree: II
- Do not use after expiration date. Do not reuse the test.
- Do not use test if foil pouch containing the test has been damaged.
- Keep out of the reach of children.
- Do not open the test foil pouch until you are ready to start the test. Once the foil pouch is opened, use within one hour.
- This IVD digital test meets the emissions and immunity requirements. The EMC countermeasures employed within the electronic instrument will provide reasonable protection against electromagnetic interference effects likely to be encountered in the home environment. Please check following preventive warnings:
 - a) Use of this instrument in a dry environment, especially if synthetic materials are present (synthetic clothing, carpets etc.) may cause damaging electrostatic discharges that may cause erroneous results.
 - b) Do not use this instrument in close proximity to sources of strong electromagnetic radiation (e.g., mobile phones), as these may interfere with the proper operation.
- This test provides only a preliminary test result, it could not replace medical advice. Always seek guidance from a qualified healthcare provider before making any health-related decisions.
- Any serious incident that has occurred in relation to the device shall be reported to the manufacturer and the competent authority.

SPECIMEN COLLECTION

The urine specimen can be collected in a clean and dry container. The first morning urine specimen is preferred; however, urine specimens collected at any time of the day may be used. Please perform the test immediately after urine collection.

DIRECTIONS FOR USE


Allow the test and specimen to reach room temperature (15–30 °C) prior to testing.

Note: Throughout testing, never hold the test with the Absorbent Tip pointing upwards.


Step 1

When you are ready, remove the test from the foil pouch and remove the cap. Hold Anti-slip Grip with the Absorbent Tip pointing downwards.

Step 2

Replace the Cap and place the test on a flat surface with the Display Screen facing upwards.

The clock symbol  will appear and blink on the Display Screen after urine is applied, indicating the test is working.

Step 3

In approximately 3 minutes, the result “Pregnant” or “Not Pregnant” will appear on the Display Screen.

READING THE RESULTS

Pregnant	shows on the display, this means that you are probably pregnant .	The result will remain on the display screen for about one hour. If there are other results shown on the display screen, please read the section of Test Errors.
Not Pregnant	shows on the display, this means that you are probably not pregnant .	

LIMITATIONS

1. The Digital Pregnancy Test is for the qualitative detection of Human Chorionic Gonadotropin in human urine, the test results are for clinical reference only and cannot be used alone as the basis for the diagnosis or exclusion of cases.
2. Drugs which contain hCG (such as Pregnyl, Profasi, APL and etc) can give a false positive result. Alcohol and painkillers that do not contain hCG should not affect the test result.
3. Very diluted urine specimens, as indicated by a low specific gravity, may not contain representative levels of hCG. If pregnancy is still suspected, a first morning urine specimen should be collected 48 hours later and tested.
4. This test may produce false positive results. A number of conditions other than pregnancy, including trophoblastic disease and certain non-trophoblastic neoplasms cause elevated levels of hCG_{1,2}. Therefore, the presence of hCG alone in urine should not be used to diagnose pregnancy unless these conditions have been ruled out.
5. This test may produce false negative results. False negative results may occur when the levels of hCG are below the sensitivity level of the test. When pregnancy is still suspected, a first morning urine specimen should be collected 48 hours later and tested. In case pregnancy is suspected and the test continues to produce negative results, see a physician for further diagnosis.
6. This test provides a presumptive diagnosis for pregnancy. A confirmed pregnancy diagnosis should only be made by a physician after all clinical and laboratory findings have been evaluated.

PRINCIPLE

The Digital Pregnancy Test is used for the qualitative detection of Human chorionic gonadotropin (hCG) in urine. During testing, the specimen reacts with the particle labeled with specific anti-hCG antibody. The mixture migrates upward on the membrane chromatographically by capillary action to react with specific capture anti-hCG antibody on the membrane and generate a colored line on test midstream. The result signal is converted to electrical signal by photodiode, and then calculated by software to display in LCD.

REAGENTS

The test contains anti-hCG antibody particles and anti-hCG antibody coated on the membrane.

STORAGE AND STABILITY

The kit can be stored at room temperature or refrigerated (2–30 °C). The test is stable through the expiration date printed on the sealed pouch. The test must remain in the sealed pouch until use. DO NOT FREEZE. Do not use beyond the expiration date.

TEST ERRORS

	An error has occurred during test. It may be that: • The absorbent tip did not keep pointing downward when applying urine. • Too little urine was applied. • Other instructions for use were not followed.
	The test did not proceed as expected. It may be that: • Urine sample was too little to start the test. • The test has reached the expiration date. • Other instructions for use were not followed.

For the both errors above, you should test with a new test following the instructions carefully.

PERFORMANCE CHARACTERISTICS

Accuracy

A clinical evaluation was conducted comparing the results obtained using the Digital Pregnancy Test to another commercially available urine hCG Rapid test. The study included 150 urine specimens, and both assays identified 50 positive and 100 negative results. The results demonstrated >99% overall accuracy of Digital Pregnancy Test when compared to the other hCG Rapid Test.

Digital Pregnancy Test	Method	Other hCG Rapid Test		Total Results
	Results	Positive	Negative	
		50	0	50
	Positive	0	100	100
	Negative	50	100	150

Sensitivity: >99.9% (92.9%~>99.9%)*
Accuracy: >99.9% (97.6%~>99.9%)*
Specificity: >99.9% (96.4%~>99.9%)*
* 95% Confidence Intervals

A clinical evaluation was conducted comparing the results obtained using the Digital Pregnancy Test to another Commercial ECLIA. The study included 100 urine specimens, and both assays identified 40 positive and 60 negative results. The results demonstrated >99% overall accuracy of Digital Pregnancy Test when compared to the other commercial ECLIA.

Digital Pregnancy Test	Method	ECLIA		Total Results
	Results	Positive	Negative	
		40	0	40
	Positive	0	60	60
	Negative	40	60	100

Sensitivity: >99.9% (91.2%~>99.9%)*
Accuracy: >99.9% (96.4%~>99.9%)*
Specificity: >99.9% (94.0%~>99.9%)*
* 95% Confidence Intervals

Sensitivity

The Digital Pregnancy Test can detect hCG levels of 25 mIU/mL or higher.

Hook Effect

The high-dose hook effect was not observed at concentrations up to 500,000 mIU/mL.

Precision – Repeatability

Precision – Repeatability was assessed using hCG concentrations of 0 mIU/mL, 12.5 mIU/mL, 15 mIU/mL, 18.75 mIU/mL, 22.5 mIU/mL, 25 mIU/mL, 50 mIU/mL, 100mIU/mL and 200 mIU/mL. The study was conducted with ten replicates per day for five consecutive days (stream testing for 5 different days, and dip testing for 5 different days) by a single operator using one lot of the Digital Pregnancy Test. No intra-lot differences were observed.

Precision– Reproducibility

Precision–reproducibility were assessed using hCG concentrations of 0 mIU/mL, 12.5 mIU/mL, 15 mIU/mL, 18.75 mIU/mL, 22.5 mIU/mL, 25 mIU/mL, 50 mIU/mL, 100mIU/mL and 200 mIU/mL. The study was conducted over five consecutive days (stream testing for 5 different days, and dip testing for 5 different days), with ten replicates per day at each of three different testing sites. Three separate lots of the Digital Pregnancy Test and three operators were used (one lot one operator per site). No difference was detected between days, sites, lots and operators.

Cross-Reactivity

The test meets the W.H.O. International Standard. No cross-reactivity was observed when adding LH (500 mIU/mL), FSH (1,000 mIU/mL), or TSH (1,000 _IU/mL) to negative (0 mIU/mL hCG) or positive (25 mIU/mL hCG) specimens.

Interfering Substance

The following potentially interfering substances were added to hCG-negative and hCG-positive specimen.

Acetaminophen	20 mg/dL	Phenothiazine	20 mg/dL	Amoxicillin	20 mg/dL
Ascorbic acid	20 mg/dL	EDTA	80 mg/dL	Uric acid	23.5 mg/dL
Atropine	20 mg/dL	Ethanol	1.0%	Estradiol valerate	20 mg/mL
Caffeine	20 mg/dL	Salicylic Acid	20 mg/dL	Ketone	20 mg/dL
Genisic acid	20 mg/dL	Bilirubin	40 mg/dL	Pregnanediol	1500 µg/dL
Glucose	2 g/dL	17β-estradiol	8 mg/mL	Phenylpropanolamine	20 mg/dL
Hemoglobin	1000 mg/dL	Ibuprofen	40 mg/dL	Benzoylcegonine	10 mg/dL
Tetracycline	20 mg/dL	Progesterone	8 mg/mL	Cannabino	10 mg/dL
Ampicillin	20 mg/dL	Vitamin B1	80 mg/dL	Ephedrine	20 mg/dL
Albumin	2000 mg/dL	Folic acid	0.03 mg/dL	Thiophene	20 mg/dL
β-hydroxybutyrate	2000 mg/dL	Aspirin	80 mg/dL		

No interference was observed from any tested substances at their respective concentrations

EXTRA INFORMATION

- 1 **How accurate is the Digital Pregnancy Test?**
A clinical evaluation of 150 urine specimens was conducted by professionals, and the results obtained through the Digital Pregnancy Test were compared with another commercial available urine hCG test. The overall accuracy is > 99%.
A clinical evaluation of 100 urine specimens was conducted by lay users, and the results obtained through the Digital Pregnancy Test were compared with another commercially ECLIA hCG test. The overall accuracy is >99%.
- 2 **How does the test work?**
It detects the pregnancy hormone hCG (human chorionic gonadotropin) in urine. When you are pregnant, your body produces the pregnancy hormone hCG (human chorionic gonadotropin), and the amount increases in the early stage of pregnancy. The Digital Pregnancy Test can detect the hormone in your urine as low as 25mIU/mL.
- 3 **Do I have to test with first morning urine?**
Although you can test at any time of the day, your first morning urine is usually the most concentrated of the day. Using morning urine to test may get more accurate result.
- 4 **What should I do if the result shows that I am “Pregnant”?**
You may be pregnant. See your doctor to confirm the result and get advice.
- 5 **What should I do if the result shows that I am “Not Pregnant”?**
You may not be pregnant or the level of pregnancy hormone may not be high enough to be detected. If you do not start your period within a week of its due date, repeat the test with a new test. If you receive the same result after repeating the test and you still do not get your period, you should see your doctor.
- 6 **There is an “?” displayed on the screen. What does it mean?**
The “?” symbol indicates an error has occurred during test, you may read the instructions carefully and test again with a new test. See Test Errors.
- 7 **I have used the test but no result has appeared on the display. What does it mean?**
If no result appears, please see Test Errors.

DISPOSING OF THE TEST

The test contains batteries, remove batteries prior to test disposal. Battery disposal should be in accordance with the appropriate recycling scheme. To remove, insert a coin into the slot at the end and twist to break open the back of the test.

Caution: Do not disassemble, recharge or dispose of batteries in fire. Dispose of the rest of the test according to the appropriate recycling scheme for electrical equipment. Do not dispose of electrical equipment in fire.

BIBLIOGRAPHY

1. Dawood MY, BB Saxena, R Landesman Human chorionic gonadotropin and its subunits in hydatidiform mole and choriocarcinoma, Obstet. Gynecol. 1977; 50(2): 172–181.
2. Braunstein GD, JL Vaitukaitis, PP Carbone, GT Ross. Ectopic production of human chorionic gonadotropin by neoplasms. Ann Intern Med 1973; 78(1): 39–45.

Index of Symbols

	Consult instructions for use or consult electronic instructions for use		Use-by date
	In vitro diagnostic medical device		Manufacturer
	Authorized representative		Humidity limitation
	Do not use if package is damaged and consult instructions for use		Temperature limit
	Dispose items according to local relevant laws regarding disposal and recycle		Catalogue number
	Importer		Do not re-use
	Contains sufficient for <n> tests		Unique device identifier
	Batch code		Recyclable materials

Importer: **Alveola Ltd.**
D-12529 Berlin, Willy-Brandt-Platz 2
H-1143 Budapest, Gizella út 28/A
hello@x-epil.com

MedNet EC-REP GmbH
Borkstrasse 10, 48163 Muenster, Germany

0123

 **Hangzhou AllTest Biotech Co., Ltd.**
#550, Yinhai Street,
Hangzhou Economic & Technological
Development Area,
Hangzhou, 310018 P.R. China
Web: www.alltests.com.cn
Email: info@alltests.com.cn

Number: **Revision date: 2026-xx-xx**